

Organogenesis

OCMP™

Osteoconductive Matrix PLUS



- Osteoconductive Matrix PLUS™ (OCMP) is a blend of ground cancellous chips and demineralized cortical powder – all derived from the same donor.
- Demineralized bone has been shown to have osteoconductive and osteoinductive potential¹
- Gamma sterilized for patient safety
- Freeze-dried for convenient ambient temperature storage with a five-year shelf life
- Available in 2cc, 5cc, and 10cc sizes for optimal efficiency

Product Code	Product Description	Dimensions
OCMP2	Osteoconductive Matrix PLUS 2cc	One (1) Package – Containing 2cc
OCMP5	Osteoconductive Matrix PLUS 5cc	One (1) Package – Containing 5cc
OCMP10	Osteoconductive Matrix PLUS 10cc	One (1) Package – Containing 10cc

References: 1. Gruskin E, Doll BA, Futrell FW, Schmitz JP, Hollinger JO. Demineralized bone matrix in bone repair: History and use. *Advanced Drug Delivery Reviews.* 2012;64(12):1063-1077.